

How to Apply ULTRABLOX® X-Ray Attenuating Cream

SINGLE GLOVE METHOD

See reverse side for Double Glove Application Method



STEP 1

Open tube and squeeze out attenuating cream into one hand.



STEP 2


Spread cream on front and back of both hands, making sure to evenly and thoroughly coat each side and between each finger. Use enough cream so that both hands are completely covered with a paper-thin layer of cream.



STEP 3

Apply surgical glove* over cream on each hand making sure to encapsulate the cream between the glove and the hand.

*For use with natural rubber latex gloves and latex-free poly-isoprene surgeon's gloves only.



REMOVAL

Wash hands using a mild abrasive soap or a scrubby bar utilizing your normal post-surgical scrub technique.



**ULTRA
BLOX**

ULTRABLOX Cream has received a 510k from the FDA. U.S. Patent No. 9,114,121.



BLOXR Solutions, LLC

P.O. Box 5148, North Branch, NJ 08876 USA

855-256-9729 | Fax 877-254-4888 | order@bloxr.com

www.bloxr.com

How to Apply ULTRABLOX® X-Ray Attenuating Cream

DOUBLE GLOVE METHOD

See reverse side for Single Glove Application Method



STEP 1

Put on the first pair of surgical gloves.



STEP 2

Open tube and squeeze out attenuating cream into gloved hands.



STEP 3

Spread cream on front and back of both gloved hands, making sure to evenly and thoroughly coat each side and between each finger. Use enough cream so that hands are completely covered with a paper-thin layer of cream.



STEP 4

Apply the second surgical glove over the cream making sure to encapsulate the cream between the two gloves.



REMOVAL

Grasp end of glove taking care to remove bottom glove, cream and top glove all at once.



**ULTRA
BLOX**
X-Ray Attenuating
Cream

ULTRABLOX Cream has received a 510k from the FDA.
U.S. Patent No. 9,114,121.



BLOXR Solutions, LLC

P.O. Box 5148, North Branch, NJ 08876 USA

855-256-9729 | Fax 877-254-4888 | order@bloxr.com

www.bloxr.com